

# SYNCHRO LINE


## BUTTERFLY-REVERSE



Butterfly-Reverse  
Art.-No. **70900** Fitness

Butterfly-Reverse  
Medical Art.-No. **7090M** Certified Medical Product

Length cm	Width cm	Height cm	Weight kg
<b>166</b>	<b>105</b>	<b>131</b>	<b>330</b>



Training of the middle fibers of the hood muscle, rear part of the deltoid muscle, large and small rhomboid muscles (M. trapezius pars transversa, M. delto-ideus pars spinalis, M. rhomboideus minor et major)

✓ Integrated

+ Optional

Adjustable starting position



Angle display



Range Limiter for End Position



Biofeedback-system



Package additional weight



### Properties

- Training in short and long lever possible
- effective stabilisation of the pectoral girdle
- specially tailored power curve
- exercise can be started conveniently by actuating the electronic boarding aid using a foot switch
- all device settings can be made fully automatic by RFID transponder according to the instructions

### Device Setting

- Precise adjustment of starting angle by using a foot switch
- max. user weight: 200kg